

Oral Fluid Drug Test Package Insert

Package insert for testing of the following drugs:

Amphetamine, Barbiturates, Benzodiazepine, Cocaine, Marijuana, Methadone, Methamphetamine, Morphine, Opiate and Propoxyphene.

For employment and insurance use.

For forensic use.

INTENDED USE & SUMMARY

The Oral Fluid Drug Test is intended for screening for the presence of drugs and their metabolites in oral fluid. For professional *in vitro* diagnostic use only.

The Oral Fluid Drug Test is a lateral flow chromatographic immunoassay for the qualitative detection of drugs and drug metabolites in oral fluid at the following cut-off concentrations:

Test	Calibrator	Cut-off (ng/mL)
Amphetamine (AMP)	d-Amphetamine	50/80
Barbiturate (BAR)	Secobarbital	50/300
Benzodiazepine (BZO)	Oxazepam	5/10/20/50
Cocaine (COC)	Benzoylcegonine	20/40
Marijuana (THC)	11-nor- Δ^9 -THC-9 COOH	12
Marijuana (THC)	Δ^9 -THC	25/50/75/100
Methadone (MTD)	Methadone	30/75
Methamphetamine (MET)	D-Methamphetamine	50
Morphine (MOP)	Morphine	15
Opiates (OPI)	Morphine	40
Propoxyphene (PPX)	Propoxyphene	50

This test will detect other related compounds, please refer to the Analytical Specificity table in this package insert.

AMP: Amphetamine is a sympathomimetic amine with therapeutic indications. The drug is often self-administered by nasal inhalation or oral ingestion.¹

BAR: Barbiturates are central nervous system depressants. They are used therapeutically as sedatives, hypnotics, and anticonvulsants. Barbiturates are almost always taken orally as capsules or tablets. The effects resemble those of intoxication with alcohol. Chronic use of barbiturates leads to tolerance and physical dependence. Short acting Barbiturates taken at 400 mg/day for 2-3 months can produce a clinically significant degree of physical dependence. Withdrawal symptoms experienced during periods of drug abstinence can be severe enough to cause death.

BZO: Benzodiazepines are central nervous system (CNS) depressants commonly prescribed for the short-term treatment of anxiety and insomnia. In general, benzodiazepines act as hypnotics in high doses, as anxiolytics in moderate doses and as sedatives in low doses. The use of benzodiazepines can result in drowsiness and confusion. Psychological and physical dependence on benzodiazepines can develop if high doses of the drug are given over a prolonged period. Benzodiazepines are taken orally or by intramuscular or intravenous injection, and are extensively oxidized in the liver to metabolites. Benzodiazepines can be detected in oral fluid after use.

COC: Cocaine is a potent central nervous system (CNS) stimulant and a local anesthetic derived from the coca plant (*Erythroxylum coca*).¹

THC: Tetrahydrocannabinol, the active ingredient in the marijuana plant (*Cannabis sativa*), is detectable in oral fluid shortly after use. The detection of the drug is thought to be primarily due to the direct exposure of the drug to the mouth (oral and smoking administrations) and the subsequent sequestering of

the drug in the buccal cavity.²

MTD: Methadone is a synthetic analgesic drug originally used for the treatment of narcotic addiction. In addition to use as a narcotic agonist, methadone is being used more frequently as a pain management agent. The psychological effects induced by using methadone are analgesia, sedation, and respiratory depression. Based on the saliva/plasma ratio calculated over salivary pH ranges of 6.4-7.6 for therapeutic or recreational doses of methadone, a cut-off <50 ng/mL is suggested. Due to this recommendation, the cut-off level of the methadone test was calibrated to 30 ng/mL.

MET: Methamphetamine is a potent stimulant chemically related to amphetamine but with greater CNS stimulation properties. The drug is often self-administered by nasal inhalation, smoking or oral ingestion.¹

OPI (MOP): The drug class opiates refers to any drug that is derived from the opium poppy, including naturally occurring compounds such as morphine and codeine and semi-synthetic drugs such as heroin. Opiates control pain by depressing the CNS, and demonstrate addictive properties when used for sustained periods of time. Opiates can be taken orally or by injection routes including intravenous, intramuscular and subcutaneous; illegal users may also take the intravenously or by nasal inhalation.³

*The window of detection varies for different opiates. Codeine can be detected within one hour and up to 7-21 hours after a single oral dose. Morphine is detectable for several days after a dose.

PPX: Propoxyphene or Dextropropoxyphene is a narcotic analgesic compound with a structural similarity to methadone. It is prescribed in the United States for the relief of moderate pain. Darvocet™, one of the most common brand names for the drug, contains 50-100 mg of propoxyphene napsylate and 325-650 mg of acetaminophen. Physiological effects of propoxyphene include respiratory depression. Propoxyphene is metabolized in the liver to yield norpropoxyphene. Norpropoxyphene has a longer half-life (30 to 36 hours) than that of propoxyphene (6 to 12 hours). Norpropoxyphene demonstrates substantially less central-nervous system depression than propoxyphene, but shows a greater local anesthetic effect.

This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) and gas chromatography/tandem mass spectrometry (GC/MS/MS) are the preferred confirmatory methods. Professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

PRINCIPLE

The Oral Fluid Drug Test is an immunoassay based on the principle of competitive binding. Drugs that may be present in the oral fluid specimen compete against their respective drug conjugate for binding sites on their specific antibody. During testing, a portion of the oral fluid specimen migrates along the test strip by capillary action. A drug, if present in the oral fluid specimen below its cut-off concentration, will not saturate the binding sites of its specific antibody. The antibody will then react with the drug-protein conjugate and a visible colored line will show up in the test line region of the specific drug strip. The presence of drug above the cut-off concentration in the oral fluid specimen will saturate all the binding sites of the antibody. Therefore, the colored line will not form in the test line region. A drug-positive oral fluid specimen will not generate a colored line in the specific test line region of the strip because of drug competition, while a drug-negative oral fluid specimen will generate a line in the test line region because of the absence of drug competition. To serve as a procedural control, a colored line will always appear at the control line region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS

The Oral Fluid Drug Test contains mouse monoclonal antibody-coupled

particles and corresponding drug-protein conjugates. A goat antibody is employed in each control line.

PRECAUTIONS

- For employment and insurance use. For forensic use.
- Do not use after the expiration date.
- The test device should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used collector and device should be discarded according to local regulations.
- Safety data sheets available for professional user upon request

STORAGE AND STABILITY

Store as packaged in the sealed pouch either at room temperature or refrigerated (2-30°C). The test device is stable through the expiration date printed on the sealed pouch. The test device must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

SPECIMEN COLLECTION AND PREPARATION

The oral fluid specimen should be collected using the collector provided with the kit. Follow the detailed Directions for Use below. No other collection devices should be used with this test. Oral fluid collected at any time of the day may be used. If specimen cannot be tested immediately, it is recommended that specimen be stored at 2-8°C or -20°C for up to 72 hours. Specimen may also be stored at room temperature for up to 48 hours. For ideal shipment conditions, transport specimen using ice packs (2-8°C).

MATERIALS

Materials Provided

- Test cups
- Saliva collectors
- Security seal labels
- Package insert

Materials Required But Not Provided

- Timer
- Gloves

DIRECTIONS FOR USE

Allow the test device, specimen, and/or controls to reach room temperature (15-30 °C) prior to testing. Instruct the donor to not place anything in the mouth including food, drink, gum, tobacco products for at least 10 minutes prior to collection.


1. Bring the pouch to room temperature before opening it. Remove the test device from the sealed pouch and use it as soon as possible.
2. Using the provided collection swab, remove the collector from the sealed pouch, have donor sweep inside of mouth (cheek, gum, tongue) several times, then hold swab in mouth until color on the saturation indicator strip appears in the indicator window of collection swab. **Important:** Do not bite, suck, or chew on the sponge.
Note: If after 7 minutes, color on the saturation indicator has not appeared in the indicator window, proceed with the test below. (See illustration 1)
3. Open the cap and place the test device on a clean and flat surface. Remove the collection sponge from the mouth and insert the sponge first into the screening device until touch the bottom of the saliva cup, pushing the cap until it locked in place of the saliva cup. **Keep upright when insert the sponge.** (See illustration 2)
4. **Test device upright on flat surface and keep upright while test is running.** Wait for the colored signal to appear in test results area. Read the results at 10 minutes.

Note: 1, Once the collection sponge locks in place, the device is airtight, tamper evident, and ready to be disposed or sent to lab for confirmation (on presumptive positive result).

2, In the case of no flowing even with enough saliva specimen, or the saliva is too thick to run, please move the device but don't tilt and keep upright back and forth on a flat and clean surface for several times. Do not tilt the device when the test is running before reading results.



INTERPRETATION OF RESULTS



INTERPRETATION OF RESULTS

(Please refer to the previous illustration)

NEGATIVE: A colored line in the control line region (C) and a colored line in the test line region (T) for a specific drug indicate a negative result. This indicates that the drug concentration in the oral fluid specimen is below the designated cut-off level for that specific drug.

***NOTE:** The shade of color in the test line region (T) may vary, but it should be considered negative whenever there is even a faint colored line.

POSITIVE: A colored line in the control line region (C) but no line in the test line region (T) for a specific drug indicates a positive result. This indicates that the drug concentration in the oral fluid specimen exceeds the designated cut-off for that specific drug.

INVALID: Control line (C) fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test device. If the problem persists, discontinue using the lot immediately and contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

- The Oral Fluid Drug Test provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) or gas chromatography/tandem mass spectrometry (GC/MS/MS) is the preferred confirmatory method.
- There is a possibility that technical or procedural errors, as well as other interfering substances in the oral fluid specimen may cause erroneous results.
- A positive test result does not indicate the concentration of drug in the specimen or the route of administration.
- A negative result may not necessarily indicate a drug-free specimen. Drug may be present in the specimen below the cut-off level of the test.
- The test does not distinguish between drugs of abuse and certain medications.

6. A positive result may be obtained from certain foods or food supplements.

PERFORMANCE CHARACTERISTICS

Analytical Sensitivity

A phosphate-buffered saline (PBS) pool was spiked with drugs to target concentrations of $\pm 50\%$ cut-off and tested with the Oral Fluid Drug Test. The results are summarized below.

Drug Conc. (Cut-off range)	AMP 50	AMP 80	BAR 50	BAR 300	BZO 5	BZO 10	BZO 20
0% Cut-off	30 0	30 0	30 0	30 0	30 0	30 0	30 0
-50% Cut-off	30 0	30 0	30 0	30 0	30 0	30 0	30 0
Cut-off	15 15	14 16	14 16	14 16	14 16	14 15	17 13
+50% Cut-off	0 30	0 30	0 30	0 30	0 30	0 30	0 30
3X Cut-off	0 30	0 30	0 30	0 30	0 30	0 30	0 30

Drug Conc. (Cut-off range)	BZO 50	COC 20	COC 40	THC 12	THC 25	THC 50	THC 75
0% Cut-off	30 0	30 0	30 0	30 0	30 0	30 0	30 0
-50% Cut-off	30 0	30 0	30 0	30 0	30 0	30 0	30 0
Cut-off	14 16	14 16	14 16	16 14	17 13	15 15	14 16
+50% Cut-off	0 30	0 30	0 30	0 30	0 30	0 30	0 30
3X Cut-off	0 30	0 30	0 30	0 30	0 30	0 30	0 30

Drug Conc. (Cut-off range)	THC 50	MTD 30	MTD 75	MET	MOP	OPI	PPX
0% Cut-off	30 0	30 0	30 0	30 0	30 0	30 0	30 0
-50% Cut-off	30 0	30 0	30 0	30 0	30 0	30 0	30 0
Cut-off	15 15	14 16	15 15	14 16	14 16	15 15	15 15
+50% Cut-off	0 30	0 30	0 30	0 30	0 30	0 30	0 30
3X Cut-off	0 30	0 30	0 30	0 30	0 30	0 30	0 30

Analytical Specificity

The following table lists the concentration of compounds (ng/mL) above which the Oral Fluid Drug Test identified positive results at 10 minutes.

AMPHETAMINE (AMP 50)		METHADONE (MTD 75)	
d-Amphetamine	50	Methadone	75
d,l-Amphetamine	125	Doxylamine	100,000
β -Phenylethylamine	4,000	Estrone-3-sulfate	100,000
Tryptamine	1,500	Phencyclidine	100,000
p-Hydroxyamphetamine	800	METHADONE (MTD 30)	
(+) 3,4-Methylenedioxyamphetamine (MDA)	150	Methadone	30
l-Amphetamine	4,000	Doxylamine	50,00
AMPHETAMINE (AMP 80)		Estrone-3-sulfate	50,00
d-Amphetamine	50	Phencyclidine	50,00
d,l-Amphetamine	125	MARIJUANA (THC 12)	
β -Phenylethylamine	4,000	11-nor- Δ^9 -THC-9 COOH	12
Tryptamine	1,500	Cannabinol	31,500
p-Hydroxyamphetamine	800	11-nor- Δ^8 -THC-9 COOH	2
(+) 3,4-Methylenedioxyamphetamine (MDA)	150	Δ^8 -THC	6,000
l-Amphetamine	4,000	Δ^9 -THC	20,000
COCAINE (COC 40)		METHAMPHETAMINE (MET)	
Benzoylcegonine	40	d-Methamphetamine	50
Cocaine	40	Fenfluramine	60,000
Cocaethylene	50	p-Hydroxymethamphetamine	400
Ecgonine	2,000	Methoxyphenamine	25,000

Ecgoninemethylester	25,000	3,4-Methylenedioxyamphetamine (MDMA)	50
N-Acetylprocainamide	18,000	l-Phenylephrine	4,000
Chlordiazepoxide	18,000	Procaine	2,000
COCAINE (COC 40)		(1R,2S)-(-) Ephedrine	400
Benzoylcegonine	20	1-Ephedrine	400
Cocaine	20	Mephentermine	800
Cocaethylene	25	(-)-Deoxyephedrine, L-Methamphetamine	3,000
Ecgonine	1,500	Ephedrine	800
Ecgoninemethylester	12,500	BENZODIAZEPINES (BZO 5)	
N-Acetylprocainamide	12,500	Oxazepam	5
Chlordiazepoxide	12,500	Alprazolam	5
MARIJUANA (THC 25)		Bromazepam	10
Δ^9 -Tetrahydrocannabinol	25	Chlordiazepoxide	10
11-nor- Δ^9 -THC-9 COOH	15	Clobazam	4
MARIJUANA (THC 50)		Clorazepate	15
Δ^9 -Tetrahydrocannabinol	50	Delorazepam	15
Δ^8 -Tetrahydrocannabinol	75	Desalkylflurazepam	15
11-nor- Δ^9 -THC-9 COOH	15	Diazepam	3
11-hydroxy- Δ^9 -THC	300	Estazolam	3
Cannabinol	2,000	Flunitrazepam	60
Cannabidiol	>10,000	α -Hydroxyalprazolam	120
MARIJUANA (THC 75)		(\pm)-Lorazepam	120
Δ^9 -Tetrahydrocannabinol	75	Midazolam	15
Δ^8 -Tetrahydrocannabinol	150	Nitrazepam	8
11-nor- Δ^9 -THC-9 COOH	15	Norchlordiazepoxide	120
11-hydroxy- Δ^9 -THC	300	Nordiazepam	15
Cannabinol	1,500	Temazepam	5
Cannabidiol	>10,000	Triazolam	15
MARIJUANA (THC 100)		BENZODIAZEPINES (BZO 10)	
Δ^9 -Tetrahydrocannabinol	100	Oxazepam	10
Δ^8 -Tetrahydrocannabinol	250	Alprazolam	6
11-nor- Δ^9 -THC-9 COOH	25	Bromazepam	12
11-hydroxy- Δ^9 -THC	500	Chlordiazepoxide	12
Cannabinol	2,500	Clobazam	6
Cannabidiol	>10,000	Clorazepate	25
BARBITURATE (BAR 50)		Delorazepam	25
Secobarbital	50	Desalkylflurazepam	25
Amobarbital	100	Diazepam	3
Alphenal	100	Estazolam	3
Aprobarbital	30	Flunitrazepam	100
Butobarbital	30	α -Hydroxyalprazolam	200
Butalbital	400	(\pm)-Lorazepam	200
Butethal	30	Midazolam	25
Cyclopentobarbital	60	Nitrazepam	12
Pentobarbital	150	Norchlordiazepoxide	200
Phenobarbital	30	Nordiazepam	25
BARBITURATE (BAR 300)		Temazepam	6
Secobarbital	300	Triazolam	25
Amobarbital	300	BENZODIAZEPINES (BZO 50)	

Alphenal	150
Aprobarbital	200
Butabarbital	75
Butalbital	2,500
Butethal	100
Cyclopentobarbital	600
Pentobarbital	300
Phenobarbital	100
OPIATE (OPI 40)	
Morphine	40
Codeine	10
Ethylmorphine	24
Hydromorphone	100
Hydrocodone	100
Levorphanol	400
Oxycodone	25,000
Morphine 3-β-d-glucuronide	50
Norcodeine	1,500
Normorphine	12,500
Nalorphine	10,000
Oxymorphone	25,000
Thebaine	1,500
Diacetylmorphine (Heroin)	50
6-Monoacetylmorphine (6-MAM)	25
Bilirubin	3,500
PROPOXYPHENE (PPX)	
Propoxyphene (PPX)	50
D-Norpropoxyphene	200

Oxazepam	50
Alprazolam	300
Bromazepam	60
Chlordiazepoxide	60
Clobazam	36
Clorazepate	125
Delorazepam	125
Desalkylflurazepam	12
Diazepam	15
Estazolam	15
Flunitrazepam	500
α-Hydroxyalprazolam	1,000
(±)-Lorazepam	1,000
Midazolam	125
Nitrazepam	60
Norchlordiazepoxide	1,000
Nordiazepam	125
Temazepam	30
Triazolam	40
MORPHINE (MOP)	
Morphine	15
Codeine	15
Ethylmorphine	15
Hydromorphone	50
Hydrocodone	50
Morphine 3-β-d-glucuronide	30
Nalorphine	300
Oxymorphone	25,000
Thebaine	5,000
Diacetylmorphine (Heroin)	15
6-Monoacetylmorphine (6-MAM)	15

Caffeine
Chloral hydrate
Chloramphenicol
Chlorothiazide
d,l-Chlorpheniramine
Chlorpromazine
Chloroquine
Cholesterol
Clonidine
Cortisone
Creatinine
Deoxycorticosterone
Dextromethorphan
Diclofenac
Dicyclomine
Diflunisal
Digoxin
Diphenhydramine
β-Estradiol
Ethyl-p-aminobenzoate
l-Epinephrine
Erythromycin
Fenopropfen
Furosemide
Gentisic acid
Hemoglobin
Hydralazine
Hydrochlorothiazide
Hydrocortisone
o-Hydroxyhippuric acid
βHydroxynorephedrine
5-Hydroxytryptamine (Serotonin)
3-Hydroxytyramine
Ibuprofen
lproniazid
(-)-Isoproterenol
Isoxsuprine
Ketoprofen
Oxymetazoline
Papaverine
Penicillin-G
Pentazocine
Perphenazine
Phenelzine
Trans-2-phenylcyclo-propylamine
Phentermine
Phenylpropanolamine
Prednisolone
Phenolbarbital
Prednisone
d,l-Propranolol
d-Pseudoephedrine
Quinacrine
Quinine
Quindine
Ranitidine
Salicylic acid
Sulfamethazine
Sulindac
Tetracycline
Tetrahydrocortisone3-acetate
Tetrahydrocortisone
3 (β-d-glucuronide)
Theophylline
Thiamine
Thioridazine
d,l-Tyrosine
Tolbutamide
Trazodone
Triamterene
Trifluoperazine
Trimethoprim
d,l-Tryptophan
Tyramine
Uric acid
Verapamil
Zomepirac

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds spiked into drug-free PBS stock. The following compounds demonstrated no false positive results on the Oral Fluid Drug Test when tested at concentrations up to 100 µg/mL.

Non Cross-Reacting Compounds

Acetaminophen	Labetalol
Acetophenetidine	Loperamide
Acetylsalicylic acid	Meprobamate
Aminopyrine	Methylphenidate
Amoxicillin	Nalidixic acid
Ampicillin	Naproxen
Amitypyline	Niacinamide
Ascorbic acid	Nifedipine
Apomorphine	Nimesulide
Aspartame	Norethindrone
Atropine	Noscapine
Benzilic acid	d,l-Octopamine
Benzoic acid	Oxalic acid
Benzphetamine	Oxolinic acid

BIBLIOGRAPHY

1. Moolchan E, *et al.* Saliva and Plasma Testing for Drugs of Abuse: Comparison of the Disposition and Pharmacological Effects of Cocaine. Addiction Research Center, IRP, NIDA, NIH, Baltimore, MD. As presented at the SOFT-TIAFT meeting October 1998.
2. Schramm W., *et al.* Drugs of Abuse in Saliva: A Review. *J Anal Tox*, 16 (1): 1-9, 1992.
3. Kim L, *et al.* Plasma and oral fluid pharmacokinetics and pharmacodynamics after oral codeine administration. *ClinChem*, 48 (9): 1486-96, 2002.
4. Kang GI and Abbott FS. Analysis of methadone and metabolites in biological fluids with gas chromatography-mass spectrometry. *J Chromatogr.* 231 (2); 311-319. Sept 1982.