

QUIKSCREEN®

ONE-STEP ON-SITE DRUG CUP
AVAILABLE WITH ALCOHOL AND ADULTERATION TEST



DRUGS OF ABUSE

Unlike other on-site tests, QuikScreen does not require any handling or manipulation to activate the test. Just dispense urine into cup and read the results.



- > Simple Procedure
- > Easy-to-interpret Color Bands
- > Built-in Test Control
- > Room Temperature Storage
- > 18-24 Month Shelf Life
- > No Instruments Needed

> Optional Adulteration Strips Built-in!

QuikScreen® is the easy-to-use urinalysis drug test that detects drugs and drug metabolites at SAMHSA cutoff levels

ORDERING INFORMATION:

No. of Parameters	Regular Cup	Adulteration Cup
17 Test Cup	60HXX	65HXX
16 Test Cup	60GXX	65GXX
15 Test Cup	60FXX	65FXX
14 Test Cup	60EXX	65EXX
13 Test Cup	60DXX	65DXX
12 Test Cup	60CXX	65CXX
11 Test Cup	60BXX	65BXX
10 Test Cup	60AXX	65AXX
9 Test Cup	609XX	659XX
8 Test Cup	608XX	658XX
7 Test Cup	607XX	657XX
6 Test Cup	606XX	656XX
5 Test Cup	605XX	655XX
4 Test Cup	604XX	654XX
3 Test Cup	603XX	653XX

"X" denotes number representation of each drug combination. Refer to catalog listing for combination and catalog number. Please contact our marketing team for the most up-to-date certification status of our drugs.

Assay Adulteration Parameters available in:

- > Alcohol
- > Creatinine
- > Glutaraldehyde
- > NPB
- > pH
- > Specific Gravity



U.S. Patent No. 6,548,019 B1

Create your own combination!

Drug Class	Abbreviations	Sensitivity
AMPHETAMINE-I	AMP	500 ng/ml
AMPHETAMINE-II	AMP (CLIA)	1000 ng/ml
BARBITURATES-I	BAR	200 ng/ml
BARBITURATES-II	BAR (CLIA)	300 ng/ml
BENZODIAZEPINES-I	BZD	200 ng/ml
BENZODIAZEPINES-II	BZD (CLIA)	300 ng/ml
BUPRENORPHINE	BUP*	10 ng/ml
COCAINE-I	COC	150 ng/ml
COCAINE-II	COC (CLIA)	300 ng/ml
COTININE	COT*	200 ng/ml
COTININE	COT*	500 ng/ml
2-ETHYLIDINE-1,5-DIMETHYL-3,3-DIPHENYLPIPERIDINE	EDDP	100 ng/ml
KETAMINE	KET*	1000 ng/ml
MARIJUANA	THC (CLIA)	50 ng/ml
METHADONE	MAD (CLIA)	300 ng/ml
METHAMPHETAMINE-I	MET	500 ng/ml
METHAMPHETAMINE-II	MET (CLIA)	1000 ng/ml
METHYLENEDIOXYMETHAMPHETAMINE-I	MDMA	500 ng/ml
METHYLENEDIOXYMETHAMPHETAMINE-II	MDMA	1000 ng/ml
OPIATES-I	OPI	300 ng/ml
OPIATES-II	OPI (CLIA)	2000 ng/ml
OXYCODONE	OXY (CLIA)	100 ng/ml
PHENCYCLIDINE	PCP (CLIA)	25 ng/ml
PROPOXYPHENE	PPX 510K	300 ng/ml
SYNTHETIC CANNABIS	K2*	50 ng/ml
TRAMADOL	TRA*	100 ng/ml
TRAMADOL	TRA*	200 ng/ml
TRICYCLIC ANTIDEPRESSANT	TCA (CLIA)	1000 ng/ml

*BUP, COT, K2, KET, TRA -- Please inquire for more information

> **ACCURACY:** FDA 510(K) cleared detection strips have the same high degree of accuracy when compared to SAMHSA certified laboratory screening methods.

> **SELF-CONTAINED:** Each cup includes a temperature strip and an identification label. Each cup is individually pouched and has an 18-24 month shelf life.

> **NO LEAKAGE:** Leak-free cups allow for positive specimen to be shipped in the same container to a laboratory for GC/MS confirmation.

> **COST:** In addition to providing immediate results, the QuikScreen® device costs less than half the price of laboratory testing.

Test Procedure

- 1 Urinate into the cup, filling the cup with urine above the minimum level line
- 2 Wait 5 minutes
- 3 Read the results
- 4 Refer to product insert for detailed instructions

QuikScreen®

One-Step Onsite Drug Cup

Confirm: A single line (top band) indicates a positive result

Negative: Two lines indicate a negative result

- Negative
- Confirm
- Retest
- Invalid