

# DRUGCHECK<sup>®</sup> SalivaScan<sup>™</sup>

**NEW groundbreaking features never before seen on the market!**

The DrugCheck<sup>®</sup> SalivaScan<sup>™</sup> oral fluid drug test includes the first sponge saturation indicator! No more throwing away devices due to inadequate saliva. Indicator strip in the collection swab changes color once the sponge contains enough saliva. The swab locks into the device for a tamper evident and airtight seal.


- **Sponge saturation indicator** (Patent Pending)
- **Market-leading drug cutoffs – No buffer required**  
**THC 50 (parent) – BUP 5 – EDDP 20**
- **100% tamper evident with saliva held in airtight chamber for confirmation**
- Advanced oral fluid collection media (sponge)
  - Sponge reduces viscosity for faster-running tests and more uniform results
  - Highest recovery of THC
  - No ascorbic acid – No additives
- Detect up to 10 drugs including alcohol
- One drug per strip for superior accuracy and sensitivity
- Easy access extraction ports for lab confirmation
- 100% made in the U.S.A.
- Crystal-clear container with magnification of the results window
- Compatible with VeriCheck<sup>™</sup> Results Manager system or photocopy/scan results
- CE marked



**FORENSIC USE ONLY**

\* Patent Pending

AVAILABLE DRUGS					
ALC ALC ALC	0.04%	COC COC COC	20 ng/mL	OPI OPI OPI	40 ng/mL
AMP AMP AMP	50 ng/mL	EDDP EDDP EDD	20 ng/mL	OXY OXY OXY	40 ng/mL
BUP BUP BUP	5 ng/mL	MET MET MET	50 ng/mL	PCP PCP PCP	10 ng/mL
BZO BZO BZO	10 ng/mL	MTD MTD MTD	30 ng/mL	THC THC THC	50 ng/mL (Parent)
Δ-9-tetrahydrocannabinol					



## MARKET-BEATING FEATURES

**1** First ever sponge saturation indicator

**2** Lowest THC cutoff: 50 ng/ml (Parent)  
 Also EDDP 20 ng/ml and BUP 5 ng/ml  
 No buffer required

**3** 100% tamper evident – Saliva held in airtight chamber for lab confirmation

Results display in same format as DrugCheck<sup>®</sup> urine devices