

T-CUP



It Is Time to Have Your Own Cup

- Female User Friendly
- No Donor Tampering Possible
- Dark Line Present As Soon As 1 Minute
- Up To 17 Drug Panels Plus Adulteration



5~14 Drug Panel Cup List

Item Code	Description
TDOA-154	5-Panel T-Cup; COC,THC,OPI,AMP,PCP
TDOA-254	5-Panel T-Cup; COC,THC,OPI,AMP,mAMP
TDOA-255	5-Panel T-Cup; COC,THC,MOP,AMP,mAMP
TDOA-454	5-Panel T-Cup; COC,THC,OPI,mAMP,PCP
TDOA-654	5-Panel T-Cup; COC,THC,OPI,mAMP,BZO
TDOA-264	6-Panel T-Cup; COC,THC,OPI,AMP,mAMP,BZO
TDOA-564	6-Panel T-cup; COC,THC,OPI,mAMP,BZO,OXY
TDOA-274	7-Panel T-Cup; COC,THC,OPI,AMP,mAMP,BZO,OXY
TDOA-474	7-Panel T-Cup; THC,COC,mAMP,OPI,PCP,AMP,MDMA
TDOA-3104	10-Panel Cup; COC,THC,OPI,AMP,mAMP,PCP,BZO,BAR,MTD,MDMA
TDOA-4104	10-Panel T-CUP; COC,THC,OPI,AMP,mAMP,PCP,BZO,BAR,MTD,OXY
TDOA-4105	10-Panel T-CUP; COC,THC,MOP,AMP,mAMP,PCP,BZO,BAR,MTD,OXY
TDOA-7104	10-Panel T-Cup; COC,THC,OPI2000,BZO,mAMP,TCA,OXY,BUP,BAR,MTD
TDOA-7105	10-Panel T-Cup; COC,THC,MOP,BZO,mAMP,TCA,OXY,BUP,BAR,MTD
TDOA-8104	10-Panel T-Cup; THC,COC,OPI,mAMP,AMP,BZO,MTD,OXY,MDMA,BUP
TDOA-3124	12-Pane T-Cup; AMP/COC/OXY/THC/PCP/MDMA/OPI/BZO/BAR/MTD/TCA/BUP
TDOA-4124	12-Panel T-Cup; COC,THC,OPI,OXY,AMP,BZO,BAR,mAMP,PPX,MTD,PCP,BUP
TDOA-4125	12-Panel T-Cup; COC,THC,MOP,OXY,AMP,BZO,BAR,mAMP,PPX,MTD,PCP,BUP
TDOA-6124	12-Pane T-Cup; AMP,COC,OXY,THC,PCP,MDMA,OPI,BZO,BAR,MTD,mAMP,BUP
TDOA-6125	12-Pane T-Cup; AMP,COC,OXY,THC,PCP,MDMA,MOP,BZO,BAR,MTD,mAMP,BUP
TDOA-1134	13-Panel T-Cup; COC,THC,OPI,BZO,mAMP,OXY,BUP,BAR,MTD,AMP,MDMA,PCP,PPX
TDOA-1144	14-Panel T-Cup; COC,THC,OPI,BZO,mAMP,TCA,OXY,BUP,BAR,MTD,AMP,MDMA,PCP,PPX
TDOA-1145	14-Panel T-Cup; COC,THC,MOP,BZO,mAMP,TCA,OXY,BUP,BAR,MTD,AMP,MDMA,PCP,PPX

Drug Test Available

510K CLEARED And CE Marked Tests:

THC, mAMP, AMP, OPI, MOP, COC, PPX, BUP, MTD, MDMA, BZO, TCA, OXY, BAR, PCP

CLIA WAIVED Test:

mAMP, OPI, MOP, COC, MTD, MDMA, TCA, PCP

(Tests in calendar: AMP, BZO, BAR, THC, BUP, OXY, PPX)

For Forensic Use only:

EDDP(300ng/mL), TRA(200ng/mL), THC(40ng/mL), AMP(500ng/mL, 300ng/mL), BZO(200ng/mL),

COC(150ng/mL, 100ng/mL), mAMP(500ng/mL), OPI(100ng/mL), MTD(200ng/mL)

Built in Adulteration parameters, offering visual check for urine integrity:

pH, Specific Gravity, Creatinine, Nitrite, Oxidant, Glutaraldehyde, Bleach, Pyridinium Chlorchromate

